


# Blom Subglottic Suctioning Cannula

## *Closing the Gap on VAP*

- Removes secretions accumulated above the cuff thereby reducing the risk of VAP
- Continuous suction at 20mmHg or Intermittent cycled suction at 100-150mmHg
- Suction line cap allows for intermittent manual suctioning when needed
- Manufactured with a suction catheter attached to the outside inner cannula which ends at the level of the fenestration
- Single Patient Use for up to 24 hrs

Valve for Intermittent Manual Subglottic Suctioning



Part No.	Product Description	Qty/Box
<b>Subglottic Suctioning Cannulas</b>		
313-50-040	Subglottic Suctioning Cannula #4 (Grey)	10
313-50-060	Subglottic Suctioning Cannula #6 (Green)	10
313-50-080	Subglottic Suctioning Cannula #8 (Blue)	10
313-50-100	Subglottic Suctioning Cannula #10 (Black)	10

**Blom**<sup>®</sup>

Complete  
line of  
Tracheostomy  
Tube products

**Pulmonary**<sup>®</sup>

*...bringing change to life*<sup>®</sup>

www.Pulmonary.com

CE  
0344

Blom is a registered trademark of Hansa Medical Products, Inc. in the United States and European Communities. Pulmonary and *bringing change to life* are registered trademarks of Pulmonary, Inc. LPV, Exhaled Volume Reservoir, EVR and SoftTouch are trademarks of Pulmonary, Inc. US Patent Pending. Patented under European Patent Number 1961441.