

Blom Fenestrated Cuffed Ordinary Tracheostomy Tubes With Extraordinary Features

- The fenestration is approximately 1mm above the cuff
- Cuff inflation prevents fenestration contact with the tracheal mucosa
- Fenestration location and rounded, smooth edges helps prevent granulation
- Used with the Standard, Fenestrated, Subglottic Suctioning and Speech Cannulas or LPV™



Cuffed Fully Deflated

Part No.	Product Description	Qty/Box
----------	---------------------	---------

Blom Tracheostomy Tube Kit (Fenestrated Cuffed Tube)

Includes Fenestrated Cuffed Tube, Standard Cannula, Obturator, Decannulation Plug and SoftTouch™

Tube Holder

313-10-040	Fenestrated Cuffed Tracheostomy Tube Kit #4	1
313-10-060	Fenestrated Cuffed Tracheostomy Tube Kit #6	1
313-10-080	Fenestrated Cuffed Tracheostomy Tube Kit #8	1
313-10-100	Fenestrated Cuffed Tracheostomy Tube Kit #10	1
313-71-000	Decannulation Plug	10
313-72-000	SoftTouch™ Tube Holder	20

Blom®

Complete
line of
Tracheostomy
Tube products

Pulmonary®

...bringing change to life®

www.Pulmonary.com

CE
0344

Blom is a registered trademark of Hansa Medical Products, Inc. in the United States and European Communities. Pulmonary and *bringing change to life* are registered trademarks of Pulmonary, Inc. LPV, Exhaled Volume Reservoir, EVR and SoftTouch are trademarks of Pulmonary, Inc. US Patent Pending. Patented under European Patent Number 1961441.