


# BLOM *UNIQUE CANNULAS*<sup>™</sup>

## MAKE OUR TRACHEOSTOMY TUBES EXTRAORDINARY



### **Disposable Subglottic Suctioning Cannula**

The Blom *Disposable* Subglottic Suctioning Inner Cannula allows continuous suctioning of secretions above the cuff for up to 24 hours.

### **Cuff-Up Speech Cannula**

Ventilator dependent patients are able to speak naturally with the Blom Speech Cannula, without cuff deflation or extensive ventilator manipulation.

### **Cosmetically Enhanced Speaking Valve**

The Blom Low Profile Valve (LPV<sup>™</sup>) can safely provide the gift of speech for tracheostomized, non-vented, spontaneously breathing patients.

*From left to right:*  
Blom Speech Cannula,  
LPV, Fenestrated Cuffed  
Tracheostomy Tube and  
Subglottic Suctioning  
Cannula

Complete  
line of  
Tracheostomy  
Tube products

For more information about **Blom Tracheostomy Tubes** and our complete line of airway products visit our website at [www.pulmodyne.com](http://www.pulmodyne.com) or call **317-246-5505**.

**Pulmodyne**<sup>®</sup>



Part No.	Product Description	Qty/Box
<b>Blom Tracheostomy Tube Kit (Fenestrated Cuffed Tube)</b>		
<i>Includes Fenestrated Cuffed Tube, Standard Cannula, Obturator, Decannulation Plug and SoftTouch™ Tube Holder</i>		
313-10-040	Fenestrated Cuffed Tracheostomy Tube Kit #4	1
313-10-060	Fenestrated Cuffed Tracheostomy Tube Kit #6	1
313-10-080	Fenestrated Cuffed Tracheostomy Tube Kit #8	1
313-10-100	Fenestrated Cuffed Tracheostomy Tube Kit #10	1
<b>Blom Tracheostomy Tube Kit (Non-Fenestrated Uncuffed Tube)</b>		
<i>Includes Non-Fenestrated Uncuffed Tube, Standard Cannula, Obturator, Decannulation Plug and SoftTouch™ Tube Holder</i>		
313-11-040	Non-Fenestrated Uncuffed Tracheostomy Tube Kit #4	1
313-11-060	Non-Fenestrated Uncuffed Tracheostomy Tube Kit #6	1
313-11-080	Non-Fenestrated Uncuffed Tracheostomy Tube Kit #8	1
313-11-100	Non-Fenestrated Uncuffed Tracheostomy Tube Kit #10	1
<b>Blom Tracheostomy Tube Kit (Fenestrated Uncuffed Tube)</b>		
<i>Includes Fenestrated Uncuffed Tube, Fenestrated Cannula, Obturator, 15mm Cap and SoftTouch™ Tube Holder</i>		
313-12-040	Fenestrated Uncuffed Tracheostomy Tube Kit #4	1
313-12-060	Fenestrated Uncuffed Tracheostomy Tube Kit #6	1
313-12-080	Fenestrated Uncuffed Tracheostomy Tube Kit #8	1
313-12-100	Fenestrated Uncuffed Tracheostomy Tube Kit #10	1
<b>Subglottic Suctioning Cannulas</b>		
313-50-040	Subglottic Suctioning Cannula #4 (Grey)	10
313-50-060	Subglottic Suctioning Cannula #6 (Green)	10
313-50-080	Subglottic Suctioning Cannula #8 (Blue)	10
313-50-100	Subglottic Suctioning Cannula #10 (Black)	10
<b>Blom Speech Cannula</b>		
313-30-040	Speech Cannula #4	10
313-30-060	Speech Cannula #6	10
313-30-080	Speech Cannula #8	10
313-30-100	Speech Cannula #10	10
<b>Blom Standard and Fenestrated Cannulas</b>		
313-20-040	Standard Cannula #4	10
313-20-060	Standard Cannula #6	10
313-20-080	Standard Cannula #8	10
313-20-100	Standard Cannula #10	10
313-22-040	Fenestrated Cannula #4	10
313-22-060	Fenestrated Cannula #6	10
313-22-080	Fenestrated Cannula #8	10
313-22-100	Fenestrated Cannula #10	10
<b>LPV (Low Profile Valve)</b>		
313-40-040	LPV #4	10
313-40-060	LPV #6	10
313-40-080	LPV #8	10
313-40-100	LPV #10	10
<b>Accessories</b>		
313-70-000	EVR™	10
313-71-000	Decannulation Plug	10
313-71-015	15mm Cap	10
313-72-000	SoftTouch™ Tube Holder	20

### Blom Tracheostomy Tubes:

- Fenestrated Cuffed
- Non-Fenestrated Uncuffed
- Fenestrated Uncuffed



### Blom Cannulas:

- Subglottic Suctioning Cannula
- Speech Cannula
- Standard Cannula
- Fenestrated Cannula
- LPV™

