

**BiTrac MaxShield™**

# **BiTrac MaxShield**

*THE NEXT GENERATION IN  
NON-INVASIVE VENTILATION*

*BILEVEL COMPATIBLE AND DISPOSABLE:*

- Prevents nasal skin breakdown
- Silicone cushion for full facial perimeter low pressure seal
- Two sizes fits most adults, simplifies mask fitting



*Complete line  
of masks for  
immediate  
non-invasive  
ventilation*

For more information about **BiTrac MaxShield**  
and our complete line of interface products visit  
our website at [www.pulmodyne.com](http://www.pulmodyne.com) or call  
**317-246-5505**.

**Pulmodyne®**

Part No.	Product Description	Size	Case/Qty.
----------	---------------------	------	-----------

**Masks**

**BiTrac MaxShield** (Bilevel and CPAP Application) w/ Anti-Asphyxia Elbow 22mm Female and MaxShield Head Strap

313-9052		Adult	10 per case
313-9054		Small Adult	10 per case

**BiTrac SE MaxShield** (Ventilator Application) w/ Standard Elbow 22mm Female and MaxShield Head Strap

313-9053		Adult	10 per case
313-9055		Small Adult	10 per case

**Circuits**

**BiTrac™ NIV Circuit** (Bilevel/CPAP Applications) 72" Smooth ID Tubing w/ swivel connector w/ exhalation port, cap plug; includes optional 24" tubing w/ "T" and proximal pressure line

313-9074			10 per case
----------	--	--	-------------

**BiTrac NIV Circuit** (Bilevel/CPAP Applications) 72" Smooth ID Tubing w/ attached filter, swivel connector w/ exhalation port, cap plug; includes optional 24" tubing w/ "T" and proximal pressure line

313-9075			10 per case
----------	--	--	-------------

**ValueTrac™ NIV Circuit** (Bilevel/CPAP Applications) 72" Corrugated Tubing and swivel connector w/ exhalation port and proximal pressure line

313-9078			10 per case
----------	--	--	-------------

Anti-Asphyxia Elbow for patient safety during Bilevel Applications (*BiTrac SE MaxShield has a Standard Elbow for ventilator application*)

BiTrac MaxShield allows unobstructed view minimizing claustrophobia

Convenient, pre-attached Head Strap

Head Strap with quick and easy attachment clips

Two exhaust ports provide continuous leak (*BiTrac SE MaxShield does not have exhaust ports*)

Manometer port for connecting proximal pressure line or for monitoring CO<sub>2</sub>



Made in USA

Single Patient Use

Rx ONLY

Latex Free



BiTrac MaxShield



BiTrac SE MaxShield



BiTrac NIV Circuit



**Pulmodyne**®

... bringing change to life®

Pulmodyne  
2055 Executive Drive  
Indianapolis, IN 46241 USA  
T: 317.246.5505  
F: 317.246.5501  
www.Pulmodyne.com

Pulmodyne and bringing change to life are registered trademarks of Pulmodyne, Inc. BiTrac, BiTrac MaxShield, Bilevel Compatible and ValueTrac are trademarks of Pulmodyne, Inc.

**CE**  
0344

8292B