

O2-RESQ™

O2-RESQ SYSTEM

*THE NEW STANDARD IN
DISPOSABLE CPAP*

SIMPLICITY AND PORTABILITY:

- For use in a wide range of respiratory distress conditions
- O₂ delivery at 30% with flow up to 140 l/m
- Easy to use **3-SET™ O2-CPAP™** Valve



*Complete line of
pre-assembled
CPAP Systems
for Emergency
Ventilation*

For more information about **O2-RESQ** and our complete line of ventilation products visit our website at www.pulmodyne.com or call **317-246-5505**.

Pulmodyne®

O2-RESQ SYSTEMS

BiTrac ED Full Face Mask with 30/22mm Elbow, OmniClip™ and Head Strap, Flow Generator and **Expandable Anti-Asphyxia Circuit** with CPAP Valve and Adapter 22mmM x 30mmM

313-7055X	BiTrac ED Mask with 3-SET O2-CPAP Valve (5.0, 7.5, 10.0cm)	Adult Large	10
313-7056X	BiTrac ED Mask with 3-SET O2-CPAP Valve (5.0, 7.5, 10.0cm)	Adult Medium	10

***Expandable Anti-Asphyxia Circuit expands to 72"**

BiTrac ED Full Face Mask with 30/22mm Elbow, OmniClip and Head Strap, Flow Generator and 72" **Corrugated Anti-Asphyxia Circuit** with CPAP Valve and Adapter 22mmM x 30mmM

313-7055	BiTrac ED Mask with 3-SET O2-CPAP Valve (5.0, 7.5, 10.0cm)	Adult Large	10
313-7056	BiTrac ED Mask with 3-SET O2-CPAP Valve (5.0, 7.5, 10.0cm)	Adult Medium	10

****Contact Pulmodyne for Additional O2-RESQ System Configurations**

ACCESSORIES

O2-CPAP™ Valves with Connectors (22mm and 30mm) 10

313-7721	20.0cm O2-CPAP Valve (Gold)	313-7717	7.5cm O2-CPAP Valve (Yellow)
313-7720	15.0cm O2-CPAP Valve (Red)	313-7716	5.0cm O2-CPAP Valve (Blue)
313-7719	12.5cm O2-CPAP Valve (Pink)	313-7715	2.5cm O2-CPAP Valve (White)
313-7718	10.0cm O2-CPAP Valve (Green)	313-7725	3-SET O2-CPAP Valve (5.0, 7.5, 10.0cm)
313-7730	5-SET™ O2-CPAP Valve (5.0, 7.5, 10.0, 12.5, 15.0cm)		

Valved "T" Adapter™ with one arm 22mmM and one arm 30mmM

313-2230			30
----------	--	--	----

O2-RESQ Filter (30mmM x 30mmF; BFE and VFE 99.9999%)

313-7043			10
----------	--	--	----

MASKS

BiTrac ED™ Full Face Mask and Head Strap

313-7031		Adult Extra Large	10
313-7030		Adult Large	10
313-7029		Adult Medium	10
313-7028		Adult Small	10

Made in USA

Single Patient Use

Rx ONLY

 Latex Free

Pulmodyne®

... bringing change to life®

Pulmodyne
2055 Executive Drive
Indianapolis, IN 46241 USA

T: 317.246.5505
F: 317.246.5501

www.Pulmodyne.com

Pulmodyne and ...bringing change to life are registered trademarks of Pulmodyne, Inc. O2-RESQ, OmniClip, 3-SET, 5-SET, O2-CPAP Valve, BiTrac ED and Valved "T" Adapter are trademarks of Pulmodyne, Inc.

CE
0344