

**Pertrach<sup>®</sup>**

# Pertrach

## *RAPID AND SIMPLE EMERGENCY AIRWAY*

*EMERGENCY CRICOTHYROTOMY OR  
EMERGENCY AND ELECTIVE TRACHEOSTOMY:*

- Patented Splitting Needle
- Adult and Pediatric Kits Available
- Adult Cuffed tube has a 5.6mm I.D.



*Innovative  
Technology  
for the Difficult  
Airway*

For more information about **Pertrach** and our complete line of airway products visit our website at [www.pulmodyne.com](http://www.pulmodyne.com) or call **317-246-5505**.

**Pulmodyne<sup>®</sup>**

**Part No.**

**Product Description**

**Case/Qty.**


**Pertrach Kits**

313-D4550	<b>Adult Pertrach, 5.6mm ID Airway, 6.8cm length, Cuffed Airway</b>	5
	Contents: Splitting needle (2) 14ga, Cuffed Trachea tube, 5.6mm ID, Dilator with flexible leader, Twill tape, 5cc Syringe, 5" Extension Tube and Scalpel	
313-D4551	<b>Same as above without Scalpel</b>	5
313-D4552	<b>Pediatric Pertrach, 3.0mm ID Airway, 3.9cm length, Uncuffed Airway</b>	5
	AGES: 6 months to 1 year (Suggested)	
	Contents: Splitting needle (2) 17ga, Uncuffed Trachea tube, 3.0mm ID, Dilator with flexible leader, Twill tape, 3cc Syringe, 5" Extension Tube and Scalpel	
313-D4553	<b>Same as above without Scalpel</b>	5
313-D4554	<b>Pediatric Pertrach, 3.5mm ID Airway, 4.0cm length, Uncuffed Airway</b>	5
	AGES: 1 to 4 years (Suggested)	
	Contents: Splitting needle (2) 17ga, Uncuffed Trachea tube, 3.5mm ID, Dilator with flexible leader, Twill tape, 3cc Syringe, 5" Extension Tube and Scalpel	
313-D4555	<b>Same as above without Scalpel</b>	5
313-D4556	<b>Pediatric Pertrach, 4.0mm ID Airway, 4.1cm length, Uncuffed Airway</b>	5
	AGES: 3 to 10 years (Suggested)	
	Contents: Splitting needle (2) 17ga, Uncuffed Trachea tube, 4.0mm ID, Dilator with flexible leader, Twill tape, 3cc Syringe, 5" Extension Tube and Scalpel	
313-D4557	<b>Same as above without Scalpel</b>	5
313-D4558	<b>Pediatric Pertrach, 5.0mm ID Airway, 4.4cm length, Uncuffed Airway</b>	5
	AGES: 10 years + (Suggested)	
	Contents: Splitting needle (2) 17ga, Uncuffed Trachea tube, 5.0mm ID, Dilator with flexible leader, Twill tape, 3cc Syringe, 5" Extension Tube and Scalpel	
313-D4559	<b>Same as above without Scalpel</b>	5


**Accessories, Educational and Training Materials**

313-A0103	5" Extension Tube	20
313-A0105	Training Dilators and Cuffed Trachea Tubes	
	Contents: 5 - Cuffed Trachea tubes, 5.6mm ID; 5 - Dilators with flexible leader; 5 - 5cc Syringes	1 Kit
313-A0107	Training Cuffed Trachea Tubes, 5.6mm ID with 5cc Syringes (10 each per kit)	1 Kit
313-A0108	Training Dilators with flexible leader (10 dilators per kit)	1 Kit
313-A0109	Metal Training Needles, 14ga (20 needles per kit)	1 Kit
313-A0110	Plastic Training Needles (20 needles per kit)	1 Kit
313-A0115	Training Scalpel Blades (10 blades per kit)	1 Kit
313-A0172	Training DVD - Pertrach insertion demonstration	1
313-A0175	Pertrach Training Block	1


**Procedural Sequence**



Simple



Rapid



Airway

Made in USA

Single Patient Use

Rx ONLY



**Pulmodyne**®

... bringing change to life®

Pulmodyne  
2055 Executive Drive  
Indianapolis, IN 46241 USA

T: 317.246.5505  
F: 317.246.5501

www.Pulmodyne.com

Pulmodyne and bringing change to life are registered trademarks of Pulmodyne, Inc. Pertrach is a registered trademark of Pulmodyne, Inc.

**CE**  
0344

6130D